

For Research Use Only

# FGF17 Polyclonal antibody

Catalog Number: 25314-1-AP

1 Publications



## Basic Information

Catalog Number:

25314-1-AP

Concentration:

850 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17807

GenBank Accession Number:

BC105131

GeneID (NCBI):

8822

UNIPROT ID:

O60258

Full Name:

fibroblast growth factor 17

Calculated MW:

216 aa, 25 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IF/ICC 1:50-1:500

## Applications

Tested Applications:

WB, IF/ICC, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Positive Controls:

WB : MCF-7 cells, PC-3 cells

IF/ICC : NIH/3T3 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yandong Gong	37172563	Cell	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 25314-1-AP (FGF17 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using FGF17 antibody (25314-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).