

For Research Use Only

FGF17 Polyclonal antibody

Catalog Number: 25314-1-AP **1 Publications**



Basic Information

Catalog Number: 25314-1-AP	GenBank Accession Number: BC105131	Purification Method: Antigen affinity purification
Size: 850 µg/ml	GeneID (NCBI): 8822	Recommended Dilutions: WB 1:500-1:1000 IF 1:10-1:100
Source: Rabbit	UNIPROT ID: O60258	
Isotype: IgG	Full Name: fibroblast growth factor 17	
Immunogen Catalog Number: AG17807	Calculated MW: 216 aa, 25 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications: IF/ICC, WB, ELISA	Positive Controls: WB : MCF-7 cells, PC-3 cells IF : NIH/3T3 cells,
Cited Applications: IF	
Species Specificity: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yandong Gong	37172563	Cell	IF


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 25314-1-AP (FGF17 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of NIH/3T3 cells using 25314-1-AP (FGF17 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).