

For Research Use Only

GNAT3 Polyclonal antibody, PBS Only

Catalog Number: 25322-1-PBS



Basic Information

Catalog Number:

25322-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18121

GenBank Accession Number:

BC147016

GeneID (NCBI):

346562

UNIPROT ID:

A8MTJ3

Full Name:

guanine nucleotide binding protein,
alpha transducing 3

Calculated MW:

354 aa, 40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

GNAT3 (Guanine nucleotide-binding protein G(t) subunit alpha-3), also known as α -gustducin, is a protein that plays a key role in taste signaling. It is part of the G-protein family and is mainly expressed in Type II taste cells in the taste buds of the tongue. GNAT3 is involved in detecting sweet, bitter, and umami tastes by activating a signaling cascade that includes PLC β 2, TRPM5, and other downstream effectors.

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse small intestine tissue slide using 25322-1-AP (GNAT3 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 25322-1-PBS in a different storage buffer formulation.