

For Research Use Only

FATP5/SLC27A5 Polyclonal antibody, PBS Only

Catalog Number: 25347-1-PBS



Basic Information

Catalog Number:

25347-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18446

GenBank Accession Number:

BC148807

GeneID (NCBI):

10998

UNIPROT ID:

Q9Y2P5

Full Name:

solute carrier family 27 (fatty acid transporter), member 5

Calculated MW:

690 aa, 75 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

SLC27A5, also known as FATP5, ACSVL6, and BAL, is an isozyme of very long-chain acyl-CoA synthetase (VLCs). It is capable of activating very long-chain fatty acids containing 24- and 26-carbons (PMID: 10479480). FATP5 mediates the import of long-chain fatty acids (LCFA) by facilitating their transport across cell membranes (PubMed:20448275). FATP5 is associated with endoplasmic reticulum but not with peroxisomes, and its primary role is in fatty acid elongation or complex lipid synthesis. FATP5 is exclusively expressed by the liver and localized to the basal plasma membrane of hepatocytes (PMID:16618416).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25347-1-AP (FATP5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25347-1-PBS in a different storage buffer formulation.



IP result of anti-FATP5 (IP:25347-1-AP, 4ug; Detection:25347-1-AP 1:5000) with rat liver tissue lysate 2695 ug. This data was developed using the same antibody clone with 25347-1-PBS in a different storage buffer formulation.