

For Research Use Only

KRT40 Polyclonal antibody

Catalog Number: 25398-1-AP



Basic Information

Catalog Number:

25398-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21889

GenBank Accession Number:

BC150191

GeneID (NCBI):

125115

UNIPROT ID:

Q6A162

Full Name:

keratin 40

Calculated MW:

431 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: A431 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 25398-1-AP (KRT40 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.