

For Research Use Only

C13orf27/TEX30 Polyclonal antibody



Catalog Number: 25408-1-AP

Basic Information

Catalog Number:

25408-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22096

GenBank Accession Number:

BC015148

GeneID (NCBI):

93081

UNIPROT ID:

Q5JUR7

Full Name:

chromosome 13 open reading frame

27

Calculated MW:

227 aa, 26 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human testis tissue, mouse testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 25408-1-AP (C13orf27/TEX30 antibody at dilution of 1:1500 incubated at room temperature for 1.5 hours.