

For Research Use Only

# PLEKHF2 Polyclonal antibody

Catalog Number: 25424-1-AP



## Basic Information

Catalog Number:

25424-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22285

GenBank Accession Number:

BC011806

GeneID (NCBI):

79666

UNIPROT ID:

Q9H8W4

Full Name:

pleckstrin homology domain  
containing, family F (with FYVE  
domain) member 2

Calculated MW:

249 aa, 28 kDa

Observed MW:

22-27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25424-1-AP (PLEKHF2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.