

For Research Use Only

NBEAL1 Polyclonal antibody, PBS Only

Catalog Number: 25431-1-PBS



Basic Information

Catalog Number:

25431-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22125

GenBank Accession Number:

BC137214

GeneID (NCBI):

65065

UNIPROT ID:

Q6ZS30

Full Name:

neurobeachin-like 1

Calculated MW:

200 aa, 307 kDa

Observed MW:

280 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Neurobeachin-like 1 (NBEAL1) is a Golgi-associated protein that plays an essential role in cholesterol metabolism regulation.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to Tris-Acetate gel system followed by western blot with 25431-1-AP (NBEAL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25431-1-PBS in a different storage buffer formulation.