

For Research Use Only

# SGLT1 Polyclonal antibody

Catalog Number: 25436-1-AP



## Basic Information

Catalog Number:

25436-1-AP

Concentration:

167 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1546

GenBank Accession Number:

BC139838

GeneID (NCBI):

6523

UNIPROT ID:

P13866

Full Name:

solute carrier family 5  
(sodium/glucose cotransporter),  
member 1

Calculated MW:

73 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

