

For Research Use Only

# BFSP2 Polyclonal antibody

Catalog Number: 25438-1-AP **1 Publications**



## Basic Information

Catalog Number:

25438-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22073

GenBank Accession Number:

BC113518

GeneID (NCBI):

8419

UNIPROT ID:

Q13515

Full Name:

beaded filament structural protein 2,  
phakinin

Calculated MW:

415 aa, 46 kDa

Observed MW:

45-49 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse eye and rat eye tissues,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jeremy A Whitson	28951044	Free Radic Biol Med	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



mouse and rat eye tissues were subjected to SDS PAGE followed by western blot with 25438-1-AP (BFSP2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.