

For Research Use Only

# SCML4 Polyclonal antibody

Catalog Number: 25439-1-AP



## Basic Information

Catalog Number:

25439-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22157

GenBank Accession Number:

BC033286

GeneID (NCBI):

256380

UNIPROT ID:

Q8N228

Full Name:

sex comb on midleg-like 4  
(Drosophila)

Calculated MW:

414 aa, 45 kDa

Observed MW:

50-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, A375 cells, HEK-293 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25439-1-AP (SCML4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.