

For Research Use Only

C2orf51 Polyclonal antibody

Catalog Number: 25464-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

25464-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21923

GenBank Accession Number:

BC029522

GeneID (NCBI):

200523

UNIPROT ID:

Q96LM6

Full Name:

chromosome 2 open reading frame 51

Calculated MW:

180 aa, 21 kDa

Observed MW:

~35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : mouse testis tissue,

Background Information

C2orf51, also named as TEX37 and TSC21, is a protein with calculated MW 21 kDa. We always detect 35 kDa in our WB detection.

Notable Publications

Author	Pubmed ID	Journal	Application
Manan Khan	29563520	Sci Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 25464-1-AP (C2orf51 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.