

For Research Use Only

# ZNF506 Polyclonal antibody

Catalog Number: 25481-1-AP



## Basic Information

Catalog Number:

25481-1-AP

Size:

950 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001099269

GeneID (NCBI):

440515

UNIPROT ID:

Q5JVG8

Full Name:

zinc finger protein 506

Calculated MW:

52 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : U-937 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 25481-1-AP (ZNF506 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.