

For Research Use Only

PUF-A, KIAA0020 Polyclonal antibody, PBS Only

Catalog Number: 25491-1-PBS

Featured Product



Basic Information

Catalog Number:

25491-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22011

GenBank Accession Number:

BC016137

GeneID (NCBI):

9933

UNIPROT ID:

Q15397

Full Name:

KIAA0020

Calculated MW:

648 aa, 74 kDa

Observed MW:

69-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of KIAA0020 antibody (25491-1-AP; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-KIAA0020 transfected A549 cells. This data was developed using the same antibody clone with 25491-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 25491-1-AP (PUF-A, KIAA0020 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25491-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human brown disease slide using 25491-1-AP (PUF-A, KIAA0020 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 25491-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human brown disease slide using 25491-1-AP (PUF-A, KIAA0020 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 25491-1-PBS in a different storage buffer formulation.