

For Research Use Only

C18orf45 Polyclonal antibody

Catalog Number: 25499-1-AP



Basic Information

Catalog Number:

25499-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21686

GenBank Accession Number:

BC114561

GeneID (NCBI):

85019

UNIPROT ID:

Q24JQ0

Full Name:

chromosome 18 open reading frame

45

Calculated MW:

296 aa, 33 kDa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 25499-1-AP (C18orf45 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.