

For Research Use Only

WDR92 Polyclonal antibody, PBS Only

Catalog Number: 25540-1-PBS



Basic Information

Catalog Number:

25540-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22181

GenBank Accession Number:

BC066657

GeneID (NCBI):

116143

UNIPROT ID:

Q96MX6

Full Name:

WD repeat domain 92

Calculated MW:

357 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




K-562 cells were subjected to SDS PAGE followed by western blot with 25540-1-AP (WDR92 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25540-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 25540-1-AP (WDR92 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 25540-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 25540-1-AP (WDR92 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 25540-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of HepG2 cells using 25540-1-AP (WDR92 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 25540-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human normal colon slide using 25540-1-AP (WDR92 antibody) at dilution of 1:1200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 25540-1-PBS in a different storage buffer formulation.