

For Research Use Only

C1orf109 Polyclonal antibody

Catalog Number: 25552-1-AP



Basic Information

Catalog Number:

25552-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22147

GenBank Accession Number:

BC018109

GeneID (NCBI):

54955

UNIPROT ID:

Q9NX04

Full Name:

chromosome 1 open reading frame

109

Calculated MW:

203 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : K-562 cells, LNCaP cells, U-251 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25552-1-AP (C1orf109 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.