

For Research Use Only

ARMC6 Polyclonal antibody

Catalog Number: 25558-1-AP



Basic Information

Catalog Number:

25558-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22278

GenBank Accession Number:

BC067125

GeneID (NCBI):

93436

UNIPROT ID:

Q6NXC6

Full Name:

armadillo repeat containing 6

Calculated MW:

501 aa, 54 kDa

Observed MW:

40-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Caco-2 cells, COLO 320 cells, HEK-293 cells, HeLa cells, Jurkat cells, U-251 cells

Background Information

ARMC6 (Armadillo repeat-containing protein 6) is a Protein Coding gene. Diseases associated with ARMC6 include Leber Congenital Amaurosis 1. The function of this gene's protein product has not been determined.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25558-1-AP (ARMC6 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 25558-1-AP (ARMC6 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.