

For Research Use Only

# DOPEY1 Polyclonal antibody, PBS Only

Catalog Number:25569-1-PBS



## Basic Information

Catalog Number:

25569-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22244

GenBank Accession Number:

BC150619

GeneID (NCBI):

23033

UNIPROT ID:

Q5JWR5

Full Name:

dopey family member 1

Calculated MW:

2456 aa, 276 kDa

Observed MW:

277 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS page followed by western blot with 25569-1-AP (DOPEY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25569-1-PBS in a different storage buffer formulation.