

For Research Use Only

ZNF323 Polyclonal antibody

Catalog Number: 25593-1-AP



Basic Information

Catalog Number:

25593-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22405

GenBank Accession Number:

BC065241

GeneID (NCBI):

64288

UNIPROT ID:

Q96LW9

Full Name:

zinc finger protein 323

Calculated MW:

406 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue, mouse lung tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 25593-1-AP (ZNF323 Antibody) at dilution of 1:600 incubated at 4 degree celsius over night.