

For Research Use Only

Znt5 Polyclonal antibody, PBS Only

Catalog Number: 25604-1-PBS

Featured Product



Basic Information

Catalog Number:

25604-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12947

GenBank Accession Number:

BC000808

GeneID (NCBI):

64924

UNIPROT ID:

Q8TAD4

Full Name:

solute carrier family 30 (zinc transporter), member 5

Calculated MW:

84 kDa

Observed MW:

84-87 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-Znt5 (IP:25604-1-AP, 4 μ g; Detection:25604-1-AP 1:500) with mouse liver tissue lysate 4000 μ g. This data was developed using the same antibody clone with 25604-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human intrahepatic cholangiocarcinoma tissue slide using 25604-1-AP (Znt5 antibody) at dilution of 1:100 (under 20 \times lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 25604-1-PBS in a different storage buffer formulation.