

For Research Use Only

EPHX3 Polyclonal antibody

Catalog Number: 25608-1-AP



Basic Information

Catalog Number:	25608-1-AP	GenBank Accession Number:	BC132960	Purification Method:	Antigen affinity purification
Size:	850 ug/ml	GeneID (NCBI):	79852	Recommended Dilutions:	IHC 1:50-1:500
Source:	Rabbit	UNIPROT ID:	Q9H6B9		IF/ICC 1:200-1:800
Isotype:	IgG	Full Name:	epoxide hydrolase 3		
Immunogen Catalog Number:	AG22277	Calculated MW:	360 aa, 41 kDa		

Applications

Tested Applications:	IHC, IF/ICC, ELISA	Positive Controls:	
Species Specificity:	human	IHC :	human breast cancer tissue,
		IF/ICC :	A549 cells,

Note-IHC: suggested antigen retrieval with **TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 25608-1-AP (EPHX3 Antibody) at dilution of 1:100 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Methanol) fixed A549 cells using EPHX3 antibody (25608-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2).