

For Research Use Only

# MOSPD2 Polyclonal antibody, PBS Only

Catalog Number: 25617-1-PBS

Featured Product



## Basic Information

Catalog Number:

25617-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22206

GenBank Accession Number:

BC030641

GeneID (NCBI):

158747

UNIPROT ID:

Q8NHP6

Full Name:

motile sperm domain containing 2

Calculated MW:

518 aa, 60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

MOSPD2 is a 518-amino acid protein containing an N-terminal CRAL/TRIO domain, a major sperm protein (MSP) domain, and a potential C-terminal transmembrane domain. MOSPD2 is a member of the vesicle-associated membrane protein-associated protein (VAP) family and one of the endoplasmic reticulum (ER) receptors that function as tether proteins to bridge ER and other organelles to build contacts (PMID: 29858488, 35389430).

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25617-1-AP (MOSPD2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25617-1-PBS in a different storage buffer formulation.



WB result of MOSPD2 antibody (25617-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-MOSPD2 transfected HeLa cells. This data was developed using the same antibody clone with 25617-1-PBS in a different storage buffer formulation.