

For Research Use Only

ZBTB22 Polyclonal antibody

Catalog Number: 25633-1-AP



Basic Information

Catalog Number:

25633-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22439

GenBank Accession Number:

BC018541

GeneID (NCBI):

9278

UNIPROT ID:

O15209

Full Name:

zinc finger and BTB domain
containing 22

Calculated MW:

634 aa, 66 kDa

Observed MW:

66-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25633-1-AP (ZBTB22 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.