

For Research Use Only

TPRG1L Polyclonal antibody

Catalog Number: 25663-1-AP



Basic Information

Catalog Number:

25663-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22585

GenBank Accession Number:

BC019034

GeneID (NCBI):

127262

UNIPROT ID:

Q5T0D9

Full Name:

tumor protein p63 regulated 1-like

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25663-1-AP (TPRG1L antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.