

For Research Use Only

LRRC29 Polyclonal antibody

Catalog Number: 25664-1-AP



Basic Information

Catalog Number:

25664-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22577

GenBank Accession Number:

BC018785

GeneID (NCBI):

26231

UNIPROT ID:

Q8WV35

Full Name:

leucine rich repeat containing 29

Calculated MW:

223 aa, 24 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, mouse pancreas tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25664-1-AP (LRRC29 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.