

For Research Use Only

PRB1 Polyclonal antibody

Catalog Number: 25675-1-AP



Basic Information

Catalog Number:

25675-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22475

GenBank Accession Number:

BC141917

GeneID (NCBI):

5542

UNIPROT ID:

P04280

Full Name:

proline-rich protein BstNI subfamily 1

Calculated MW:

392 aa, 39 kDa

Observed MW:

39-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25675-1-AP (PRB1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.