

For Research Use Only

IL-1RL2 Polyclonal antibody

Catalog Number: 25677-1-AP **1 Publications**



Basic Information

Catalog Number: 25677-1-AP	GenBank Accession Number: BC107064	Purification Method: Antigen affinity purification
Size: 240 µg/ml	GeneID (NCBI): 8808	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q9HB29	
Isotype: IgG	Full Name: interleukin 1 receptor-like 2	
Immunogen Catalog Number: AG21974	Calculated MW: 575 aa, 65 kDa	
	Observed MW: 65 kDa	

Applications

Tested Applications: FC, WB, ELISA	Positive Controls: WB : K-562 cells,
Cited Applications: FC	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
WonKyung J Cho	36740182	Am J Pathol	FC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 25677-1-AP (IL1RL2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



1×10^6 K-562 cells were stained with 0.2ug IL1RL2 antibody (25677-1-AP, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.