

For Research Use Only

# IL-1RL2 Polyclonal antibody

Catalog Number: 25677-1-AP

1 Publications



## Basic Information

Catalog Number:

25677-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21974

GenBank Accession Number:

BC107064

GeneID (NCBI):

8808

UNIPROT ID:

Q9HB29

Full Name:

interleukin 1 receptor-like 2

Calculated MW:

575 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : K-562 cells,

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
WonKyung J Cho	36740182	Am J Pathol	FC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 25677-1-AP (IL1RL2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.