

For Research Use Only

WDR38 Polyclonal antibody

Catalog Number: 25683-1-AP



Basic Information

Catalog Number: 25683-1-AP	GenBank Accession Number: BC150645	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 401551	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q5JTN6	
Isotype: IgG	Full Name: WD repeat domain 38	
Immunogen Catalog Number: AG22444	Calculated MW: 314 aa, 34 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls:
Species Specificity: human	IHC : human breast cancer tissue, human skeletal muscle tissue
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 25683-1-AP (WDR38 Antibody) at dilution of 1:50 (under 40x lens).