

For Research Use Only

TTC28 Polyclonal antibody, PBS Only

Catalog Number: 25701-1-PBS

Featured Product



Basic Information

Catalog Number:

25701-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22465

GenBank Accession Number:

BC016465

GeneID (NCBI):

23331

UNIPROT ID:

Q96AY4

Full Name:

tetratricopeptide repeat domain 28

Calculated MW:

2481 aa, 271 kDa

Observed MW:

240-271 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

TTC28 is a large protein required during the cell cycle for condensation of spindle midzone microtubules, formation of the midbody, and completion of cytokinesis. RT-PCR ELISA detected highest TTC28 expression in adult ovary and lung, followed by kidney and liver. Lower expression was detected in all other adult and fetal tissues examined except testis, which showed no expression.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25701-1-AP (TTC28 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25701-1-PBS in a different storage buffer formulation.



WB result of TTC28 antibody (25701-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-TTC28 transfected HeLa cells. This data was developed using the same antibody clone with 25701-1-PBS in a different storage buffer formulation.



IP result of anti-TTC28 (IP:25701-1-AP, 4ug; Detection:25701-1-AP 1:2500) with HEK-293 cells lysate 1480 ug. This data was developed using the same antibody clone with 25701-1-PBS in a different storage buffer formulation.