

For Research Use Only

# C17orf64 Polyclonal antibody

Catalog Number: 25706-1-AP



## Basic Information

<b>Catalog Number:</b> 25706-1-AP	<b>GenBank Accession Number:</b> BC048806	<b>Purification Method:</b> Antigen affinity purification
<b>Concentration:</b> 400 ug/ml	<b>GeneID (NCBI):</b> 124773	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q86WR6	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 17 open reading frame 64	
<b>Immunogen Catalog Number:</b> AG22769	<b>Calculated MW:</b> 236 aa, 27 kDa	
	<b>Observed MW:</b> 27-30 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse liver tissue,
<b>Species Specificity:</b> human, mouse	

## Background Information

CHCT1 (CHD1 helical C-terminal domain-containing protein 1), also called C17orf64, may play a role in the regulation of apoptosis.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25706-1-AP (C17orf64 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.