

For Research Use Only

TTC7B Polyclonal antibody

Catalog Number: 25713-1-AP



Basic Information

Catalog Number:

25713-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22388

GenBank Accession Number:

BC153122

GeneID (NCBI):

145567

UNIPROT ID:

Q86TV6

Full Name:

tetratricopeptide repeat domain 7B

Calculated MW:

843 aa, 94 kDa

Observed MW:

50-53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25713-1-AP (TTC7B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.