

For Research Use Only

TTC7B Polyclonal antibody

Catalog Number: 25713-1-AP



Basic Information

Catalog Number: 25713-1-AP	GenBank Accession Number: BC153122	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 145567	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q86TV6	
Isotype: IgG	Full Name: tetra-tricopeptide repeat domain 7B	
Immunogen Catalog Number: AG22388	Calculated MW: 843 aa, 94 kDa	
	Observed MW: 50-53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25713-1-AP (TTC7B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.