

For Research Use Only

# AGPAT3 Polyclonal antibody

Catalog Number: 25723-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

25723-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22855

GenBank Accession Number:

BC011971

GeneID (NCBI):

56894

UNIPROT ID:

Q9NRZ7

Full Name:

1-acylglycerol-3-phosphate O-acyltransferase 3

Calculated MW:

376 aa, 43 kDa

Observed MW:

36-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yanmei Peng	38383581	Exp Mol Med	WB
Shiwu Zhang	37899701	J Cachexia Sarcopenia Muscle	WB
Sasha M Milbeck	32602905	Invest Ophthalmol Vis Sci	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 25723-1-AP (AGPAT3 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.