

For Research Use Only

# PRIC285 Polyclonal antibody

Catalog Number: 25736-1-AP



## Basic Information

Catalog Number:

25736-1-AP

Concentration:

800 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22719

GenBank Accession Number:

BC150283

GeneID (NCBI):

85441

UNIPROT ID:

Q9BYK8

Full Name:

peroxisomal proliferator-activated  
receptor A interacting complex 285

Calculated MW:

2649 aa, 295 kDa

Observed MW:

294 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total  
protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: PC-3 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



IP result of anti-PRIC285 (IP:25736-1-AP, 4ug;  
Detection:25736-1-AP 1:800) with PC-3 cells lysate  
1200 ug.