

For Research Use Only

OATP5A1/SLCO5A1 Polyclonal antibody

Catalog Number: 25744-1-AP



Basic Information

Catalog Number:

25744-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22620

GenBank Accession Number:

BC137424

GeneID (NCBI):

81796

UNIPROT ID:

Q9H2Y9

Full Name:

solute carrier organic anion
transporter family, member 5A1

Observed MW:

110-12 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : DU 145 cells, mouse brain tissue, PC-3 cells, rat
brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25744-1-AP (SLC05A1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 25744-1-AP (SLC05A1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.