

For Research Use Only

TRNP1 Polyclonal antibody

Catalog Number: 25748-1-AP



Basic Information

Catalog Number:	25748-1-AP	GenBank Accession Number:	BC069216	Purification Method:	Antigen affinity purification
Size:		GeneID (NCBI):	388610	Recommended Dilutions:	IHC 1:20-1:200
Source:	Rabbit	UNIPROT ID:	Q6NT89		
Isotype:	IgG	Full Name:	TMF1-regulated nuclear protein 1		
Immunogen Catalog Number:	AG22876	Calculated MW:	227 aa, 23 kDa		

Applications

Tested Applications:	IHC, ELISA	Positive Controls:	IHC : mouse brain tissue,
Species Specificity:	human, mouse		
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0			

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 25748-1-AP (TRNP1 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 25748-1-AP (TRNP1 Antibody) at dilution of 1:50 (under 40x lens).