

For Research Use Only

TR2IT1 Polyclonal antibody

Catalog Number: 25753-1-AP



Basic Information

Catalog Number:

25753-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22682

GenBank Accession Number:

BC130548

GeneID (NCBI):

645840

UNIPROT ID:

Q6F5E7

Full Name:

thioredoxin reductase 2 intronic transcript 1

Observed MW:

18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 25753-1-AP (TR2IT1 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.