

For Research Use Only

# TR2IT1 Polyclonal antibody

Catalog Number: 25753-1-AP



## Basic Information

<b>Catalog Number:</b> 25753-1-AP	<b>GenBank Accession Number:</b> BC130548	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 400 µg/ml	<b>GeneID (NCBI):</b> 645840	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q6F5E7	
<b>Isotype:</b> IgG	<b>Full Name:</b> thioredoxin reductase 2 intronic transcript 1	
<b>Immunogen Catalog Number:</b> AG22682	<b>Observed MW:</b> 18 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : HeLa cells,
<b>Species Specificity:</b> human	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 25753-1-AP (TR2IT1 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.