

For Research Use Only

CCDC102A Polyclonal antibody

Catalog Number: 25764-1-AP



Basic Information

Catalog Number:

25764-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22753

GenBank Accession Number:

BC008285

GeneID (NCBI):

92922

UNIPROT ID:

Q96A19

Full Name:

coiled-coil domain containing 102A

Calculated MW:

550 aa, 63 kDa

Observed MW:

63-65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 25764-1-AP (CCDC102A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.