

For Research Use Only

TOM1 Polyclonal antibody

Catalog Number: 25813-1-AP



Basic Information

Catalog Number:

25813-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22931

GenBank Accession Number:

BC005144

GeneID (NCBI):

10043

Full Name:

target of myb1 (chicken)

Calculated MW:

492 aa, 54 kDa

Observed MW:

50-54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, mouse heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25813-1-AP (TOM1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.