

For Research Use Only

ZMAT2 Polyclonal antibody

Catalog Number: 25819-1-AP



Basic Information

Catalog Number:

25819-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22909

GenBank Accession Number:

BC056668

GeneID (NCBI):

153527

UNIPROT ID:

Q96NCO

Full Name:

zinc finger, matrin type 2

Calculated MW:

199 aa, 24 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : Jurkat cells, HeLa cells, SH-SY5Y cells, C6 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25819-1-AP (ZMAT2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 25819-1-AP (ZMAT2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.