

For Research Use Only

ACD Polyclonal antibody

Catalog Number: 25849-1-AP

2 Publications



Basic Information

Catalog Number:

25849-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23021

GenBank Accession Number:

BC016904

GeneID (NCBI):

65057

UNIPROT ID:

Q96AP0

Full Name:

adrenocortical dysplasia homolog
(mouse)

Calculated MW:

58 kDa

Observed MW:

57-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : Jurkat cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ji Hoon Lee	34321211	Sci Adv	IP
Hongmiao Hu	37624885	Sci Adv	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 25849-1-AP (ACD Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.