

For Research Use Only

# SOCS1 Polyclonal antibody

Catalog Number: 25852-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 25852-1-AP	<b>GenBank Accession Number:</b> BC153088	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 210 µg/ml	<b>GeneID (NCBI):</b> 8651	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> O15524	
<b>Isotype:</b> IgG	<b>Full Name:</b> suppressor of cytokine signaling 1	
<b>Immunogen Catalog Number:</b> AG23119	<b>Calculated MW:</b> 211 aa, 24 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB, IHC

**Species Specificity:**  
human

**Cited Species:**  
human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yanting Shi	34262937	Front Mol Biosci	WB, IHC
Li Jiang	32782536	Oncol Lett	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data