

For Research Use Only

SOCS1 Polyclonal antibody

Catalog Number: 25852-1-AP

3 Publications



Basic Information

Catalog Number:

25852-1-AP

Concentration:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23119

GenBank Accession Number:

BC153088

GeneID (NCBI):

8651

UNIPROT ID:

O15524

Full Name:

suppressor of cytokine signaling 1

Calculated MW:

211 aa, 24 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yanting Shi	34262937	Front Mol Biosci	WB, IHC
Li Jiang	32782536	Oncol Lett	WB
Chuan-Xia Zhang	30816302	Cell Death Differ	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

