

For Research Use Only

MTFMT Polyclonal antibody

Catalog Number: 25854-1-AP



Basic Information

Catalog Number:

25854-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23064

GenBank Accession Number:

BC033687

GeneID (NCBI):

123263

UNIPROT ID:

Q96DP5

Full Name:

mitochondrial methionyl-tRNA
formyltransferase

Calculated MW:

389 aa, 44 kDa

Observed MW:

43-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta lysates were subjected to SDS PAGE followed by western blot with 25854-1-AP (MTFMT antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.