

For Research Use Only

PHF23 Polyclonal antibody, PBS Only

Catalog Number: 25862-1-PBS

Featured Product



Basic Information

Catalog Number:

25862-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23044

GenBank Accession Number:

BC002509

GeneID (NCBI):

79142

UNIPROT ID:

Q9BUL5

Full Name:

PHD finger protein 23

Calculated MW:

403 aa, 44 kDa

Observed MW:

~65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

PHD finger protein 23 (PHF23), also known as hJUNE-1b, belongs to the PHF23 family. As a putative epigenetic reader, the PHD domain of PHF23 specifically binds to tri-methylated lysine 4 of histone H3 (H3K4me3). PHF23 is widely expressed in human tissues and various cell lines and exerts a negative effect on autophagy by promoting the degradation of E3 ligase LRSAM1 (Leucine-rich repeat and sterile alpha motif containing 1) (PMID: 32978226, 32380078). The observed molecular weight of PHF23 is ~65 kDa, which may be caused by post-translational modification. The specificity of this antibody was demonstrated by siRNA-mediated knockdown of target protein.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25862-1-AP (PHF23 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25862-1-PBS in a different storage buffer formulation.



WB result of PHF23 antibody (25862-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PHF23 transfected HeLa cells. This data was developed using the same antibody clone with 25862-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using PHF23 antibody (25862-1-AP) at dilution of 1:200 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-Phalloidin (red). This data was developed using the same antibody clone with 25862-1-PBS in a different storage buffer formulation.