

For Research Use Only

Collagen Type I Polyclonal antibody

Catalog Number: 25870-1-AP

3 Publications



Basic Information

Catalog Number:

25870-1-AP

Concentration:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22942

GenBank Accession Number:

BC036531

GeneID (NCBI):

1277

UNIPROT ID:

P02452

Full Name:

collagen, type I, alpha 1

Calculated MW:

1464 aa, 139 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

mouse, rat, sheep

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Li-Xia Qiu	39913466	PLoS One	WB
Jiancheng Lin	39722825	Front Med (Lausanne)	WB, IHC, IF
Fengqin Wang	36654751	ACS Pharmacol Transl Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

