

For Research Use Only

PIK3C2G Polyclonal antibody

Catalog Number: 25904-1-AP



Basic Information

Catalog Number:

25904-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23177

GenBank Accession Number:

BC130277

GeneID (NCBI):

5288

UNIPROT ID:

O75747

Full Name:

phosphoinositide-3-kinase, class 2,
gamma polypeptide

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB: A549 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25904-1-AP (PIK3C2G antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.